

## Why Buy MTI?

### 1. Heightened Care Experience

Provide patients and caregivers with the care experience they'll love by using MTI's innovative ADA compliant chairs featuring easy patient access, enhanced caregiver ergonomics, and smooth and ultra-quiet movements which set the patient at ease throughout the care experience.

### 2. Improve Patient Outcomes

Better outcomes are more attainable when care is provided using the world's most advanced equipment from MTI which help patients feel safe, comfortable and at ease.

### 3. Patient and Caregiver Safety

Low entry/exit heights on MTI's ADA compliant chairs ease access for all patients and empower patients with disabilities to transfer themselves while reducing the risk of injury and fatigue to healthcare workers.

### 4. Low Total Cost of Ownership

Reduce downtime and operating costs with durable MTI products that are designed for quick and easy service and repairs to optimize efficiency.

### 5. Highest Quality on the Market

MTI manufactures only the highest quality products, providing long-lasting equipment rated to ensure reliable care for even the heaviest patients and high volume use.

### 6. Affordable Innovation

MTI has launched 11 next-generation, technically innovative chairs in the last five years with features like patented SmartTech™ and ComfortSync™ technologies, which are often lower cost than outdated competitive products.

### 7. Designed for the Specialist

Based on the feedback from specialists, MTI has designed a complete offering of exam/procedure room equipment with the features and functions specialists really need.

## 528/529 Podiatry/Wound Care Chairs



Seat Height Range 18.75" - 39" (529), 21" fixed (528)  
 Seat Width 24"  
 Patient Capacity 725 lbs  
 Programmable 4 Positions  
 Base Options Standard or swivel  
 Controls 4-function (529), 3-function (528), programmable foot control, optional hand control  
 Legrest Design Fixed angle w/11" powered extension  
 Standard Armrests Removable slide-back patient arms



### Standard Features Include:

- ADA Compliant\* (529)
- Power lift (not on 528)
- Power back
- Power tilt
- 11" independent power foot extension
- Aseptic stainless steel debris tray
- SmartTech™ and ComfortSync™ technology
- Removable upholstery
- Soft Start/Stop
- Programmable home
- ADA Compliant home
- G2 Accessories & Colors


\* Only when purchased with the available transfer support accessory

725 lb Patient Capacity

## 529W Bariatric Podiatry/Wound Care Chair



Seat Height Range 18.75" - 39"  
 Seat Width 30"  
 Patient Capacity 850 lbs  
 Programmable 4 Positions  
 Base Options Standard  
 Controls 4-function programmable foot control, optional hand control  
 Legrest Design Fixed angle w/11" powered extension  
 Standard Armrests Removable slide-back patient arms



### Standard Features Include:

- ADA Compliant\*
- Power lift
- Power back
- Power tilt
- 11" independent power foot extension
- Aseptic stainless steel debris tray
- SmartTech™ and ComfortSync™ technology
- Removable upholstery
- Soft Start/Stop
- Programmable home
- ADA Compliant home
- G2 Accessories & Colors


\* Only when purchased with the available transfer support accessory

850 lb Patient Capacity

## 550 Podiatry/Wound Care Chair



Seat Height Range 18.75" - 39"  
 Seat Width 25"  
 Patient Capacity 725 lbs  
 Programmable 4 Positions  
 Base Options Standard, swivel or mobile  
 Controls 4-function programmable foot control, optional hand control  
 Legrest Design Independent footrest & powered leg extension  
 Standard Armrests Removable slide-back patient arms



### Standard Features Include:

- ADA Compliant\*
- Removable upholstery
- Power lift • Foot Accessory Splines
- Power back • Soft Start/Stop
- Power tilt • Programmable home
- Power foot w/7" independent power extension • G2 Accessories & Colors
- SmartTech™ and ComfortSync™ technology
- Rectangular non-articulating headrest

\* Only when purchased with the available transfer support accessory

725 lb Patient Capacity

## 840 Stretcher Bed



Seat Height Range 18.75" - 37"  
 Seat Width 25"  
 Patient Capacity 725 lbs  
 Programmable 4 Positions  
 Base Options Mobile  
 Controls Includes sealed membrane hand control  
 Legrest Design Independent footrest & powered leg extension  
 Standard Armrests Multiple optional arm configurations



### Standard Features Include:

- ADA Compliant\*
- Deep Oval or Ophthalmic
- Power lift articulating headrest
- Power back • 2 Lithium-ion batteries & a • Power tilt desktop charger
- Power foot w/7" independent power extension • Soft Start/Stop
- Battery or AC powered • Programmable home
- Lockable steering casters • ADA Compliant home
- SmartTech™ and • G2 Accessories & Colors
- ComfortSync™ technology

\* Only when purchased with the available transfer support accessory

725 lb Patient Capacity

## 721 Surgery Chair



Seat Height Range 19" - 35"  
 Seat Width 21.5"  
 Patient Capacity 650 lbs  
 Programmable 4 Positions  
 Base Options Standard, swivel or mobile  
 Controls Choice of either hand or foot control  
 Legrest Design Comfortable fixed angle  
 Standard Armrests Multiple optional arm configurations



### Standard Features Include:

- Power lift base with 16" of vertical travel and a 19" entry substrates height • Multiple IV/Surgery and
- Power back patient arms available
- Power tilt • G1 Accessories & Colors

650 lb Patient Capacity



Specialized Equipment for Specialty Healthcare  
 Medical Technology Industries



MTI helps specialty physicians improve patient outcomes, increase caregiver safety and efficiency, and reduce operating costs through intuitively-designed, ADA-compliant and feature-rich specialized medical and dental exam room equipment.



## MTI Products

### Cabinets

- Types of Cabinets: 9
- Countertop Height: 33.5 - 41.5"
- Width: 25-35.5"
- Drawers: 2-6
- Finish: laminate or stainless steel • Counter: laminate or solid surface



### Instrument Tables

- Locking casters
- Counter Finish: laminate or solid surface
- Counter Size: 23.75" x 39" or 24.75" x 40"
- 17.5" vertical adjustment
- Optional pull-out tray



### Stools

- Types of Stools: 12
- Controls: hand or foot
- Height: 18.25-28.75"
- Seat Diameter: 16.5"
- Backrest Options: none, ergo, curved body or twin arms
- Weight Capacity: 225-350 lbs.



### Side Chairs

- Types of Chairs: 2
- Seat Height: 17.75"
- Seat Width: 18.4"
- Backrest: ergonomic
- Weight Capacity: 425 lbs. • Arms: with or without
- Adjustable leveling feet



### Lights


- Types of Lights: 3
- Model: Clarus™ or Lumerus™ LED • Lux: 21,500 - 33,400
- Life: 40,000 - 70,000 hrs
- CRI: ≥ 75
- Style: caster base, floor stand, or chair-mounted



## 829/830 Procedure Chairs



Seat Height Range 18.75" - 39"  
 Seat Width 22" (829) or 25" (830)  
 Patient Capacity 725 lbs  
 Programmable 4 Positions  
 Base Options Standard, swivel or mobile  
 Controls Includes both hand and foot controls  
 Legrest Design Independent footrest & powered leg extension  
 Standard Armrests Multiple optional arm configurations



### Standard Features Include:

- ADA Compliant\* • Removable upholstery
- Power lift • Soft Start/Stop
- Power back • Programmable home
- Power tilt • ADA Compliant home
- Power foot w/7" independent • G2 Accessories & Colors power extension
- SmartTech™ and ComfortSync™ technology
- Oval articulating headrest


\* Only when purchased with the available transfer support accessory

725 lb Patient Capacity

## 463/464 Exam/Procedure Chairs



Seat Height Range 18.75" - 37"  
 Seat Width 21"  
 Patient Capacity 725 lbs  
 Programmable 4 Positions  
 Base Options Swivel  
 Controls EasyView™ chair back mounted membrane switches  
 Legrest Design Leg rest rises as back is reclined  
 Standard Armrests Removable slide-back patient arms



### Standard Features Include:

- ADA Compliant\*
- Power lift
- Power back & ComfortSync™ coordinated foot
- Power tilt (only on 464)
- SmartTech™ and ComfortSync™ technology
- Oval articulating headrest • Removable upholstery
- Soft Start/Stop
- Programmable home
- ADA Compliant home
- G2 Accessories & Colors


\* Only when purchased with the available transfer support accessory

725 lb Patient Capacity

## 461 Exam Chair



Seat Height Range 21" - 38"  
 Seat Width 21"  
 Patient Capacity 550 lbs  
 Programmable 2 Positions  
 Base Options Swivel top  
 Controls EasyView™ chair back mounted membrane switches  
 Legrest Design Leg rest rises as back is reclined  
 Standard Armrests Flip-up floating arms



### Standard Features Include:

- Power lift
- Power back and coordinated foot
- 21" wide seat
- Coordinated back & leg rest
- Oval articulating headrest
- Removable upholstery
- Flip-up floating patient arms • Soft Start/Stop
- 2 dynamic programmable positions
- Home button
- G2 upholstery colors

550 lb Patient Capacity

## 403 Exam Chair



Seat Height Range 21" - 38"  
 Seat Width 21"  
 Patient Capacity 550 lbs  
 Programmable 2 Positions  
 Base Options Swivel top  
 Controls EasyView™ chair back mounted membrane switches  
 Legrest Design Fixed looped footrest  
 Standard Armrests Flip-up floating arms



### Standard Features Include:

- Power lift
- Power back
- 21" wide seat
- Fixed leg position
- Oval non-articulating headrest
- Removable upholstery
- Flip-up floating patient arms • Soft Start/Stop
- 2 dynamic programmable positions
- Home button
- G2 upholstery colors

550 lb Patient Capacity

## ENT Cabinets



Cabinet Models 5  
 Height Range 33" - 41.5"  
 Depth Options 18" or 21"  
 Drawers 5 - 6  
 Doors 2-3  
 Width 25" or 35.5"  
 Pressure Systems Internal or central pressure  
 Vacuum Systems Standard or high volume internal or central vacuum  
 Finish Options Laminate or stainless steel

### Standard Features Include:

- Industrial-rated pressure and vacuum system
- Control panel with easy to use single switch vacuum/pressure operation and auxiliary switch control
- Rounded and milled side panels and edges
- Concealed heavy-duty door hinges
- Dual accessory electrical outlets located inside cabinet
- Lower cabinet bumper for safety and durability
- Latex safe (latex-free)
- Total component modularity with plug and play components
- Waste drop hole on wider models to aid in the safe sterile disposal of waste
- Configurable countertops
- Dual rechargeable instrument wells or transformer with corded handles
- Choice of 2 instrument heads • Right or left hand hose dex terit y

## Contact Us Today



**Geceg Surgical**  
 Empowering Healthcare, Enhancing Lives



**Authorized Distributor for sales, service, and support, Contact US :+1 (317) 804-1810**